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YK240 Corticosterone EIA

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FOR RESEARCH LABORATORY USE ONLY

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VIII. 文献

1. Vázquez-Palacios G, Retana-Marquez S et al: **Further definition of the effect of Corticosterone on the sleep-wake pattern in the male rat.** *Pharmacol Biochem Behav* 70:305-310, 2001
2. Goymann W, Möstl E and Gwinner E: **Corticosterone Metabolites can be Measured Noninvasively in Excreta of European Stonechats (*Saxicola torquata rubicola*).** *The Auk* 119:1167-1173, 2002
3. Hupe JM, James AC et al: **Cortical feedback improves discrimination between figure and background by V1, V2 and V3 neurons.** *Nature* 394:784-787, 1998
4. Kitaysky AS, Kitaiskaiya EV et al: **Dietary restriction causes chronic elevation of corticosterone and enhances stress response in red-legged kittiwake chicks.** *J Comp Physiol B: Regulatory, Integrative and Comparative Physiology* 171:701-709, 2001
5. Thelling O, Noel G et al: **Stress hormone secretion and gut signal transducer (STAT) proteins after burn injury in rats.** *Shock* 16:393-397, 2001
6. National Committee for Clinical Laboratory Standards Evaluation Protocols, SC1 (1989), Vallanov, PA: NCCLS
7. Takaya J, Iharada A, et al: **Upregulation of Hepatic 11β-Hydroxysteroid Dehydrogenase-1 Expression in Calcium-Deficient Rats.** *Ann Nutr Metab.* 59:73-78, 2011
8. Satoh S, Suzuki H et al: **Casted-immobilization downregulates glucocorticoid receptor expression in rat slow-twitch soleus muscle.** *Life Sciences.* 89:962-967, 2011
9. Ekuni D, Tomofuji K et al: **Occlusal disharmony increases amyloid-β in the rat hippocampus,** *NeuroMolecular Medicine,* 13:197-203, 2011
10. Irie K, Ekuni D et al: **Occlusal disharmony induces BDNF level in rat submandibular gland.** *Archives of Oral Biology.* 56:35-40, 2011

11. Yamawaki Y, Fuchikami M et al: **Antidepressant-like effect of sodium butyrate (HDAC inhibitor) and its molecular mechanism of action in the rat hippocampus.** *World Journal of Biological Psychiatry.* 2011 Aug 3. [Epub ahead of print]
12. Sato S, Suzuki H et al: **Casted-immobilization downregulates glucocorticoid receptor expression in rat slow-twitch soleus muscle.** *Life Sciences.* 89:962–967; 2011
13. Takei Y, Bartolo RC et al: **Water deprivation induces appetite and alters metabolic strategy in *Notomys alexis*: unique mechanisms for water production in the desert.** *Proc Bio Sci,* 2012 Mar 7. [Epub ahead of print]
14. Hashizume Y, Shirato K et al: **Diallyl Disulfide Reduced Dose-Dependently the Number of Lymphocyte Subsets and Monocytes in Rats.** *J Nutri Sci and Vitamin,* 58:292-296; 2012
15. Tanokashira D, Morita T et al: **Glucocorticoid Suppresses Dendritic Spine Development Mediated by Down-Regulation of Caldesmon Expression.** *J Neurosci,* 32:14583-14591; 2012
16. Nawata Y, Kitaiichi K, and Yamamoto T: **Increases of CRF in the amygdala are responsible for reinstatement of methamphetamine-seeking behavior induced by footshock.** *Pharmacol Biochem and Behav.* 101:297-302, 2012
17. Takaya J, Iharada A, et al: **Down-regulation of hepatic phosphoenolpyruvate carboxykinase expression in magnesium-deficient rats.** *Megnesium Research.* 25:131-139, 2012
18. Tomita M, Okuyama T et al: **Cardiotoxicity of methamphetamine under stress conditions: Comparison of single dose and long-term use.** *Mol Med Reports.* 7:1786-1790; 2013
19. Tanaka Y, Nakano J et al: **Hindlimb Suspension Does Not Influence Mechanical Sensitivity, Epidermal Thickness, and Peripheral Nerve Density in the Glabrous Skin of the Rat Hind Paw.** *Physiol Res,* 62:119-123; 2013
20. Takaya J, Iharada A, et al: **A calcium-deficient diet in pregnant, nursing rats induces hypomethylation of specific cytosines in the 11 β -hydroxysteroid dehydrogenase-1 promoter in pup liver.** *Nutrition Research.* 33:961-970;2013
21. Ekuni D, Endo Y et al: **Effects of apoE Deficiency and Occlusal Disharmony on Amyloid-Beta Production and Spatial Memory in Rats.** *Plos/one.* Published: September 16, 2013; DOI: 10.1371/journal.pone.0074966
22. Norii-Hayashi N, Sasagawa T et al: **Developmental Changes in Desensitisation of c-Fos Expression Induced by Repeated Maternal Separation in Pre-Weaned Mice.** *J Neuroendocrinol.* 25:158-167; 2013
23. Utsugi C, Miyazono S et al: **Hard-Diet Feeding Recovers Neurogenesis in the Subventricular Zone and Olfactory Functions of Mice Impaired by Soft-Diet Feeding.** *Plos/one.* Published:

May 9, 2014; DOI: 10.1371/journal.pone.0097309

24. Iida S, Hara T et al: **Memory-related gene expression profile of the male rat hippocampus induced by teeth extraction and occlusal support recovery.** *Arch Oral Biol.* 59:133-141, 2014
25. Goriki A, Hatanaka F et al: **A Novel Protein, CHRONO, Functions as a Core Component of the Mammalian Circadian Clock.** *Plos/Biology.* Published: April 15, 2014; DOI:10.1371/journal.pbio.1001839
26. Matsusue Y, Horii-Hayashi N et al: **Distribution of Corticosteroid Receptors in Mature Oligodendrocytes and Oligodendrocyte Progenitors of the Adult Mouse Brain.** *J Histochem Cytochem.* 62:211-226; 2014
27. Matsuda K, Horikawa Y et al: **The Adenosine Receptor Agonist 5'-N-Ethylcarboxamide-Adenosine Increases Glucose 6-Phosphatase Expression and Gluconeogenesis.** *Pharmacol & Pharmacy,* 5:19-23; 2014
28. Takaya J, Yamanouchi S et al: **A Calcium-Deficient Diet in Rat Dams during Gestation and Nursing Affects Hepatic 11 β -hydroxysteroid dehydrogenase-1 Expression in the Offspring.** *Plos/one.* Published: January 10, 2014; DOI: 10.1371/journal.pone.0084125

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