

Anti Carassius RFamide Serum

Cat. No. Y470

Lot No. 530081218

Description: This antiserum was raised in a rabbit by immunization with synthetic peptide of Carassius RFamide. The product vial contains 50 μ L of the titled antiserum obtained by lyophilizing its 0.001 M phosphate buffer (pH 7.0, 0.5mL) solution. It can be used for immunoassay, immunohistochemistry or other immunoreactive applications with Carassius RFamide.

Immunogen: Synthetic Carassius RFamide, conjugate free

Host: Rabbit

Amino Acid Sequence of RFamide (Carassius auratus)^{1,2)}:

SPEIDPFWYV GRGVRPIGRF-NH₂

Product Form: Lyophilized unpurified serum

Size: 50 μ L

Reconstitution: Reconstitute the product with 0.5mL of 0.01M PBS (pH 7.0) to make a 10 fold diluted stock solution. If it is stored in a refrigerator, add moderate antiseptic to the solution (e.g. NaN₃ 0.1%).

Storage: The product will be stable for over one year if it be stored at -20°C to -80°C until opened. Upon reconstitution, the antiserum solution must be stored at 2°C to 8°C and used within one month. Repeated freezing-thawing should be avoided.

Suggested Working Dilution Range: 1: 1,000~5,000 for immunohistochemistry (frozen or paraffin sections). 1:15,000~ for enzyme immunoassay. Optimal dilution should be determined by each laboratory for each application.


Specificity (based on non-competitive EIA): Carassius RFamide 100%, rat PrRP (12-31)-NH₂ (rat PrRP20) 100%, human PrRP (12-31)-NH₂ (human PrRP20) 27.3%, rat PrRP (3-31)-NH₂ 41.2%, human PrRP (3-31)-NH₂ 22.5%.

Tissue tested (by immunohistochemistry): with frog and carp brain giving clear positive immunostainings.

REFERENCES:

- 1) M. Fujimoto, K. Takeshita et al., Isolation and characterization of a novel bioactive peptide, Carassius RFamide (C-RFa), from the brain of the Japanese crucian carp. *Biochem. Biophys. Res. Commun.* 342(2):436-40, 1998
- 2) H. Satake, H. Minakata et al., Characterization of a cDNA encoding a precursor of Carassius RFamide, structurally related to a mammalian prolactin-releasing peptide. *FEBS Lett.* 446(2-3):247-50, 1999

FOR RESEARCH LABORATORY USE ONLY
DO NOT USE ORGANIC SOLVENTS FOR DISSOLVING ANTISERUM


Yanaihara Institute Inc.

2480-1 Awakura, Fujinomiya, Shizuoka, 418-0011 Japan
FAX: +81-544-22-2770 TEL: +81-544-22-2771

For more information, please visit our internet website: <http://www.yanaihara.co.jp> or contact us to ask@yanaihara.co.jp