

**BACKGROUND**

Rat, mouse and rabbit tumor necrosis factor -  $\alpha$  (TNF-  $\alpha$ ) are a polypeptide of 156 amino acid residues; The apparent molecular masses of TNF-  $\alpha$  under denaturing conditions are 17KDa.<sup>23)</sup> Current evidence suggests that the membrane anchored form of TNF-  $\alpha$  that exists on the surface of macrophages and/or monocytes in addition to serving as a reservoir for release of soluble TNF-  $\alpha$  has cytotoxic activity.<sup>49)</sup>

**IMMUNOGEN:** Recombinant TNF-  $\alpha$  (Rat)

**ANTIBODY TYPE:** IgG1/k

**for ELISA:** Approximately 0.4-1.2ng/well of recombinant rat TNF-  $\alpha$  can be detected using an antiserum concentration of x1000.

**for WESTERN BLOTTING:**

An antibody concentration of dilution ratio x1000 will allow visualization of 30-10ng/lane of recombinant rat TNF-  $\alpha$  under reducing condition.

**SPECIFICITY:** TNF-  $\alpha$  (Rat) 100%, TNF-  $\alpha$  (Human) <0.05%  
IL-1  $\alpha$  (Rat) < 0.05%, IL-1  $\alpha$  (Human) < 0.05%  
IL-1  $\beta$  (Rat) < 0.05%, IL-1  $\beta$  (Human) < 0.05%  
Other cytokines < 0.05%

**RELATED ANTISERA:** Rabbit Anti TNF-  $\alpha$  (Rat) Serum YC 030  
Goat Anti TNF-  $\alpha$  (Rat) Serum YC 031

**RELATED PEPTIDES:** TNF-  $\alpha$ , TNF-  $\beta$ , Other Cytokines

**STORAGE:** Keep frozen below -20°C  
Avoid repeated freezing-thawing

**REFERENCES:**

- 1) J.Vilcek and T.H. Lee, J.Biol.Chem.,266: 7313,1991
- 2) B.B. Aggarwal, B.Moffat, R.N.Harkins. J.Biol.Chem., 259: 686,1984
- 3) B.B. Aggarwal, W.J.Kohr, P.E.Hass, B.Moffat, S.A.Spencer, W.J.Henzel, T.S.Bringman, G.E.Nedwin, D.V. Goeddel, R.N.Harkins. J.Biol.Chem. 260: 2345,1985
- 4) M.Kriegler, C.Perez, K.DeFay, I.Albert, S.D.Lu. Cell. 53: 45,1988
- 5) B.Luettig, T.Decker, M.L.Lohmann-Matthes. J.Immunol. 143: 4034,1989

**FOR RESEARCH LABORATORY USE ONLY**

DO NOT USE ORGANIC SOLVENTS FOR DISSOLVING ANTISERUM